

# 30 sec Quick Wipes ANIOS

Large spectre  
d'EFFICACITÉ  
Broad spectrum  
EFFICACY



- Disinfectant wipes for non-immersible medical devices
- No rinse
- Ethanol based
- Wide spectrum in short time
- Aldehyde and colorant free formula
- Pleasant fragrance

## INDICATIONS

ANIOS QUICK WIPES are impregnated with a fast-acting disinfectant solution for non-immersible medical devices, previously cleaned and alcohol-resistant (stethoscopes, cables and connectors, pressure sensors, blood sugar testers, etc.).

## COMPOSITION

Impregnation solution : Ethanol (55%), Quaternary ammonium propionate, fragrance.

## CHARACTERISTICS

- Composition: cellulose, viscose, acrylic binder (spunlace 23g/m<sup>2</sup>)
- Wipe dimensions: 130 x 190 mm
- Impregnation solution: 1.9 ml of ANIOSPRAY QUICK per wipe.

## PHYSICOCHEMICAL DATA

The ANIOS QUICK WIPES impregnation solution is:

- Clear and colourless
- Fragrance : lemon tea
- Density at +20°C: ≈ 0.9
- pH: 5.75 +/- 0.25

# ANIOS Quick Wipes

## INSTRUCTIONS FOR USE

- Apply the wipe over the previously cleaned medical device.
- Respect the contact time indicated for the antimicrobial activity required.
- No need to rinse unless the medical device comes into contact with the skin or mucosa.
- Well close the packaging after each opening.

Do not reuse the wipe. Disinfection would no longer be guaranteed.

Comments: Before application, check the compatibility of medical devices treated with alcohol. Electrical devices for disinfection must be unplugged.



## MICROBIOLOGICAL PROPERTIES

The impregnation solution is:

Active against	Standards	Contact time
<b>Bacteria</b>	EN 1040, EN 13727, EN 14561, EN 13697. MRSA according to EN 14561	<b>30 seconds</b>
<b>Mycobacteria</b>	Mycobacterium terrae : EN 14348, EN 14563	<b>30 seconds</b>
	Mycobacterium avium : EN 14348, EN 14563	<b>5 minutes</b>
<b>Yeasts</b>	Candida albicans : EN 1275, EN 13624, EN 14562, EN 13697.	<b>30 seconds</b>
<b>Moulds</b>	Tricophyton mentagrophytes according to EN 14562	<b>30 seconds</b>
	Aspergillus niger : EN 1275, EN 13624, EN 14562, EN 13697.	<b>5 minutes</b>
<b>Virus</b>	Adenovirus, HIV-1, PRV (surrogate of HBV), BVDV (surrogate of HCV), Vaccinia virus* Herpes-virus, Coronavirus, Rotavirus, Influenza virus A (H <sub>1</sub> N <sub>1</sub> )	<b>30 seconds</b>
	Calicivirus Felin MNV according to EN 14476	<b>5 minutes</b>

\*RKI: The Robert Koch Institute recommends the use of disinfectants with the claim "begrenzt viruzid wirksan", i.e. limited viricidal activity, to interrupt the transmission chain. According to this institute, limited viricidal activity is activity of the vaccine and BVDV on the virus.

## PRECAUTIONS FOR USE

Dangerous – respect the precautions for use (Drawn up according to the European rules in force regarding the classification and labelling of chemical products).  
Storage : between +5°C and +25°C.  
Class IIb medical device (Directive 93/42/EEC as amended).

## PACKAGING

Box of 12 dispensers, 120 wipes each,  
Wipe size : 130 x 190 mm.....Ref : 2333.421

