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*Escherichia coli* WDCM 00013 (negativo)

**IMPORTANT: NEW INCUBATION TIME**

Oxidase Sticks are subject to the following modifications, due to an improvement of the raw materials:

- Uninoculated sticks present a slight pinkish appearance.
- Reaction time has been slowed (60 seconds).
- The test presents less false positives.

**Oxidase Sticks for microbiology**

For the rapid detection of cytochrome oxidase enzyme


**Product No. 416444.2326**

**Contents of the Kit**

Oxidase Sticks for microbiology -- 50 sticks

**Procedure**

1. Take the strip and wait until it reaches room temperature.
2. Select a well developed colony in the culture. Sample using the inoculation loop.
3. Using the inoculation loop, spread the colony around the test area.
4. Evaluate the result immediately. If the test area turns blue before the first **60 seconds** it means that cytochrome-oxidase test is positive. After this time all the strains, even the negative ones, cause a colour change to bluish grey.

**Precautions**

- Remove only as many strips as are required for testing.
- **Reseal the container** immediately after removing the test strips.
- Do not touch the test pads.
- It is recommended that **Nutritive Agar or Soya Tryptone Agar** are used as these do not contain colourings, indicators or inhibitors.
- **Inoculation loops containing steel, nichrome or iron** must not be used as these may result in false positives due to oxidation. A platinum or plastic loop, or wooden applicator stick is recommended.
- **Colonies** which have been extracted from a medium with a **pH of less than 5.5** must not be examined as these may result in false negatives. In the event of this, it is recommended that the colony is subjected to an intermediate culture where it is impossible for the pH level to drop below the recommended value.

**Storage and disposal**

- Protect the strips from sunlight and humidity.
- Store the package in the refrigerator (at 2-8 °C).
- Check the expiry date before the analysis.
- After use, the strips must be treated as contaminated material and must be disposed of according to good microbiological practice.
- Strips may be disposed of using combustion, autoclaving or immersing in a disinfecting solution for a minimum of 6 hours.

**Quality control**

Bacteria used as a control: *Pseudomonas aeruginosa* WDCM 00025 (positive)  
*Escherichia coli* WDCM 00013 (negative).