




IRIS *Salmonella*®

DETECTION OF ALL *SALMONELLAE* IN FOOD AND FEED PRODUCTS AND ALSO IN ENVIRONMENTAL SAMPLES



PERFORMANCE

Detection of all motile and non-motile *Salmonellae* including atypical species and serovars

EASY

Excellent contrast between clear agar and **PINK-MAGENTA** *Salmonella* colonies and the secondary flora which appears in **BLUE** to **DARK-BLUE**

RAPID

Negative screening within 37 hours and confirmation test in two minutes with a latex kit : **CONFIRM'***Salmonella*

FLEXIBLE AND ECONOMIC

Easily adaptable to all laboratory organization :
1 supplemented enrichment broth + 1 plate

RELIABLE

Validated method by AFNOR Certification according to the NF EN ISO 16140-2 and including samples greater than 25 g



IRIS *Salmonella*®

IRIS *Salmonella*® method allows the rapid detection of all *Salmonellae* in food and feed products and also in environmental samples.



Scan me



Validated by AFNOR Certification under the reference BKR 23/07-10/11

SAMPLE ≤ 25 g

SAMPLE > 25 g¹

(x) g of sample in 9 (x) mL of
Salmonella Enrichment medium²

+

Add selective **IRIS *Salmonella*® supplement**



⊕ Incubation
16 to 24 h
41.5 ± 1.0°C

⊕ Incubation
18 to 24 h
41.5 ± 1.0°C

D + 1

D + 1

Plate 10 µL on **IRIS *Salmonella*® Agar**

⊕ Incubation
24 ± 3 h
37 ± 1°C

⊕ Incubation
24 ± 3 h
37 ± 1°C

D + 2

D + 2

Reading and confirmation

Non characteristic colonies :
ABSENCE of *Salmonella*

Characteristic magenta colonies:
PRESUMPTIVE PRESENCE of *Salmonella*



Latex agglutination test :
CONFIRM' *Salmonella*

To know

¹IRIS *Salmonella* method is also validated for the detection of *Salmonellae* in samples from 50 g to 375 g in milk powder (including infantile formula with and without probiotics) and from 50 g to 125 g for flours and croquettes in feedstuffs.

² For the samples greater than 25g, preheat the *Salmonella* Enrichment medium at 41.5°C.

Salmonella Enrichment and **IRIS *Salmonella*® Agar** can be stored up to 3 days at 2-8°C before plating or reading (except the enrichment for the sample greater than 25 g in feedstuffs).

Salmonella Enrichment can be used for the preparation of dilutions.

IRIS *Salmonella*® Agar may be used as the second medium according to NF EN ISO 6579.

Please refer to the technical data sheet for more information.

Ordering information

Salmonella Enrichment

BK194HA/GC – 500 g bottle/5 kg drum
BM13608 – 10 vials of 225 mL
BM13708 – 3 flexible bags 3 L
BM14408 – 2 flexible bags 5 L

Salmonella Enrichment + Tween®80

BM16308 – 3 flexible bags 3 L
BM19808 – 2 flexible bags 5 L
BM22808 – 10 vials of 225 mL

Salmonella Enrichment double-strength buffered

BM20008 – 2 flexible bags 5 L
BM20108 – 10 vials of 225 mL
BK225HA/GC – 500 g bottle/5 kg drum

IRIS *Salmonella*® supplement

BS07808 – 10 vials of 50 mL
BS07708 – 120 tablets qsp 225 mL
BS09308 – 120 tablets qsp 90 mL

IRIS *Salmonella*® concentrated supplement

BS10208 – 10 vials of 90 mL **NEW**

IRIS *Salmonella*® Agar

BM16008 – 20 Petri dishes (Ø90mm)
BM16108 – 120 Petri dishes (Ø90mm)
BK212HA – 500 g bottle

CONFIRM' *Salmonella*

BT01108 – Latex agglutination kit (50 tests)